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EN Directions for Use – Itch Relief MosQuit® – Electronic device against insect bites

How to use your MosQuit® • For the optimum effectiveness of your MosQuit®, you should use the device as early as possible after you have been bitten. Later application, however, can also effectively combat the negative consequences of an insect bite. The heating surface of your insect bite healer is located under the nose of the mouse. Place the heating surface exactly over the puncture. Keep your insect bite healer slightly pressed on during the time of application. The MosQuit® is an innovative high-tech device that has proven itself for years on the market. The function of the device is constantly being optimized. The choice of the right temperature range, and the determination of the individual application time can have a very positive impact on the effect of the MosQuits®. There are **two temperature ranges** which can be selected for the application. The **lower temperature range** (from approx. 45°) is for children or adults with sensitive skin. The **higher temperature range** (from approx. 50° C) is particularly suitable for adults with a normal skin. If the lower temperature should be selected, the blue button above the LED has to be pressed once. Then the LED flashes in green colour until the application temperature is achieved. To select the higher temperature, the blue button above the LED has to be pressed twice. Then the LED flashes in red colour until the higher application temperature is achieved. After reaching the selected application temperature the LED glows constantly and the unit gives acoustic signals at intervals of about 1 second. Now the heating surface of the MosQuit® should be pushed lightly on the puncture. The duration of application can be determined with the help of the acoustic signals of the device. The selected application temperature should develop its positive effect on the puncture for approx. 3 – 5 seconds. After this time the MosQuit® should be removed from the puncture. The application duration should be 5 seconds for adults with normal skin and 3 seconds for children or adults with sensitive skin. The treatment can be aborted at any time by lifting the unit from the injection site (for example, when the heating is considered to be too hot or the device stop to give beeps). Since the optimal application time can be different for each user, the MosQuit® does not work with a fixed duration of application. The user himself can determine the optimal duration time by turning off the equipment at any time. During the application of the MosQuit® in the higher temperature range the unit can be switched off at once by pressing the blue button for one time. In the lower temperature range (green LED) a pressing of the blue button for one time leads to a change to the higher temperature range or the unit can be switched off at once by pressing the blue button twice.

Frequency and time of application • If possible, use the MosQuit® as early as possible after you have been bitten to reach an optimal effect. With this unit you can do several applications directly one after the other for example at different parts of the skin.

Side effects • For proper usage no side effects are known. There will be no risk of tissue damages such as burns, necrosis (cellular changes) or systemic effects when the device is used on sound and normal skin. However, an unpleasant heat feeling can be the result. For this reason, you should really apply your insect bite healer to insect bites only. Children and adults who have a sensitive skin should select at the beginning always the short heating period and the lower temperature level.

Contraindications • Your MosQuit® is no children's toy and not suitable for application to children below three years. Due to the heat it produces, do not use the MosQuit® within the delicate areas of your eyes and lips and on mucous membranes. Also, refrain from treating pre-affected skin or open wounds. The MosQuit® is no medical emergency measure in case of well-known or acute allergy to insect poisons of horseflies, bees or wasps and can only avoid or considerably alleviate itching, swelling or the hazard of infection. In case of any unusual complaints after an insect bite, consult a doctor immediately. Nor does this MosQuit® substitute any prophylactic malaria treatment for stays abroad. Diabetics should not use the MosQuit®, because they may have a reduced pain sensitivity and would be able to inflict burns on themselves. It is not allowed to use the MosQuit® by in taking of alcohol, drugs or pain reducing pharmaceuticals. The MosQuit® should only be used on one person because of hygienically reasons.

Changing the battery • Before a complain the battery has to be proven and has to be changed if it is necessary. The MosQuit® is powered by a conventional 9 V compound battery. To put in/change the battery please open the battery compartment lid on the back of the device. Put in a new battery or replace the exhausted battery and slide back the lid over the battery compartment until it snaps in.

Cleaning and Unit Disposal • Always remove the battery from the device before you clean it. The MosQuit® can be cleaned possibly with soap sud lightly moistened cloth. Do not use any chemical cleaners or abrasive agents for cleaning. Make sure, that no water comes into the device. Otherwise it may cause malfunctions or a damage of the electronics. If the water penetrates please only use the device again once it is completely dry. Please dispose the unit according to the "Electro- and Electronic Altgeräte Verordnung 2002/96/EC-WEEE". Empty batteries should not be thrown into the domestic waste. You are legal obliged to dispose the batteries correctly. If there are further inquiries please contact the communal public authority which is responsible for the disposal.

Specifications:

Dimensions: 9 cm x 5 cm x 4 cm • Heating temperature: 45 °C - 55 °C • Weight: 30 g without 9 V compound battery • Power supply: 9 V compound battery

The battery is not included in delivery.

Abbreviation:

 Read the user manual •  Do not dispose via domestic waste •  Technical changes without prior notice • Revision level 11091611

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The Declaration of Conformity is in the download center on our website www.davita.de and can be downloaded

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