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**Light therapy box**

**Appareil de lumineothérapie**

**Lichttherapie apparaat**

## **vitality Collagen**

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**User manual**

**Mode d'emploi**

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Deutsch.....	2
English.....	6
Français.....	10
Nederlands.....	14

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Dear customers,

We are happy that you have chosen a product from our company and congratulate you on this choice. You have purchased a product that can contribute to help improve your quality of life for years.

## **Effect of light therapy**

With the beginning of the age of 25 the production of collagen decreases. This leads to an accelerated skin aging, wrinkles, cellulite, loosing skin and uneven skin tone. Collagen is one of the major proteins of the body and collagen fibers are the most common components of connective tissue.

The **collagen light therapy units** use special wavelengths which achieve exactly the skin layers in which the connective tissue fibres' are produced. The collagen production can be stimulated and the collagen in the connective tissues can be increased. So the application of the light therapy can lead to an anti-aging effect with a visible rejuvenation of the skin. Wrinkles and fine lines can be refilled. Pores can be refined, the metabolism can be activated and the whole skin can be improved. The connective tissue can be tightened considerably and the elasticity of the skin can be increased significantly.

## **How to use the DAVITA® light therapy unit**

The light therapy unit can be used for the Anti-Aging-application to improve the collagen level, to reduce wrinkles, fine lines, blemishes and cellulite. The skin elasticity can be increased.

## **Some basic information before using the light therapy unit**

Before use, kindly check that the instrument is in functional order and that there is no sign of damage. Phototherapy instruments must only be operated in dry and closed rooms that are free of any explosion risks. The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide. The instruments must be standing free and should not be exposed to sun.

Before using the unit the parts of the skin should be cleaned with a pH-neutral soap and then should be carefully dried. If the units are used, there should be no rests of cosmetic/sun cream at the skin of the user.

## **How to carry out light therapy**

Place the unit on a flat surface such as a table and connect the plug with a wall socket.

The Collagen-Light therapy unit vitality Collagen has on the back side a fixture for folding, so that the device has a save stand. Sit down in front of the instrument and keep a distance of approx. 25 cm. The application duration is in this distance approx. 30 minutes.

Move the tilting switch to „On“ and the device begins to glow. After the application time the switch has to be switched off manually. During the application you can read, watch TV or carry out any other activities. A constant staring into the light source must be avoided.

Since this device emits a harmless light, it is not necessary to set up for this application a protection goggle. If you not like to look in a coloured light, it is possible to wear a protection goggle during the application.

If you want to cover the phototherapy instrument after use allow the instrument to cool down at least 15 minutes.

**Important:** In order to be protected against UV-rays, the light box must only be operated with its original protective filter in place.

## **Frequency and time of application, side effects**

The application should be carried out daily for best results or minimum 3 – 5 times per week for 30 minutes until the skin appearance has improved. Subsequently, as required, as soon as the symptoms appear again.

By using phototherapy the following side effects have been observed: Headaches, nausea, irritability and burning of the eyes can occur rarely for a short period. An absolute contraindication that means circumstances that do not allow the use of light therapy units is not known so far.

Please take care in case of all types of eye diseases and the eye specialist should be consulted before. In case of patients being under treatment with lithium or medicine for example Fluoxetine, Imipramine et al., which increase the sensibility of the eyes on light a supervision of the light therapy by an eye specialist could be considered.

These possible side effects / contraindications listed above can be avoided, if a protection goggle is worn during the application. If there is as side effect an itchy skin rash further treatments should be made only in consultation with the dermatologist.

## Fastening of the unit on the wall

The vitality Collagen can be mounted on the wall.

On the back side of the top of this unit there is an opening in the middle below the air vents (**See illustration**).

At this opening the unit can be fixed on the wall with a screw. Depending on the type of wall, an appropriate screw must be used, may be with pegs. E.g. a round head wood screw which is screwed into an expansion anchor. The screw and the anchor are not included in the scope of delivery.

## Operate the change of the lamp

Before replacing the lamp, the plug has to be disconnected from the socket. In the lower area on the back of the unit there are 2 openings. With a screwdriver 2 screws have to be loosen (**See figure 1**). Then the front part of the housing can be removed. The two tube transparent retaining clips can also be removed with a screwdriver. Then the tube can carefully be pulled out of the housing. (**See Figure 2 & 3**). After inserting the new lamp, the assembly has to be done in accordance to the instructions above but in reverse order.

Fig. 1

Fig. 2

Fig. 3

## Cleaning, maintenance and recycling

The instrument does not need any special maintenance and requires no calibration. The manufacturer does not require such preventive inspections by other persons. Qualified companies / individuals can make repairs or the change of the lamp by the device.

The mains have to be disconnected, before opening the device. To clean the outer casing parts of the device, we recommend a mild cleaning agent (as used for glass).

Packaging: DAVITA® takes part in the system „Grüner Punkt“ („Green Point“). All packaging parts are made out of recyclable materials and must be recycled in accordance of the guideless concerning each material (eco-friendly).

The casing and the front of the cover are made of plastic material. Fluorescent lamps as well as instruments marked with the symbol of a crossed out dustbin may not be disposed via the regular garbage.

The electrical and electronic components must be recycled separately. Please get in touch with your local authority to find the right way to of disposal. If electrical appliance is disposed of in landfills or dumps, hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.

## Technical data

<b>Instrument Type:</b>	<b>Light therapy box vitality Collagen</b>
<b>Article number:</b>	20450
<b>Housing:</b>	High quality plastic housing
<b>Housing measures:</b>	H: 27 cm, W: 21 cm, D: 8 cm
<b>Weight:</b>	0,8 kg
<b>Lamps:</b>	1 x 45 Watts
<b>Light spectrum:</b>	633 nm
<b>Lighting power, distance and therapy length:</b>	5.000 Lux at 25 cm distance for 30 minutes
<b>Average lifetime of the lamp via the specification of the producer:</b>	Approx. 1.000 hours
<b>Electric ballast unit:</b>	Instant starting, flickerfree light, no stroboscopic effects. Buzz free operation. Automatic switch off at the end of lamp live time.
<b>Light emission:</b>	the selection of lighting source and light filter reduces the energy of light emitted into the infra-red and UV range to a minimum level, so it is harmless to human skin and eyes.
<b>Storage:</b>	-20 degrees C up to 50 degrees C and a relative humidity of max. 90%
<b>Operation in the altitude of max.:</b>	2.000 Meter
<b>Operation temperature and humidity:</b>	10° C to 25° C; 30% to 90% relative humidity
<b>Mains requirements:</b>	220 - 240 V~ , 50 Hz

## Signs and symbols

The following symbols appear in this manual, on the device or on its accessories:

Caution: Consult accompanying documents

CE Mark

Class II equipment

**IP20**

Protected against solid foreign objects, not protected against water

Manufacturer

**SN**

Serial number

**DISPOSAL:** Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.

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Technical changes can be applied at any time without prior notice. Rev. Nr.: 11071511

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DAVITA® Medizinische Produkte GmbH & Co. KG • Im Hammereisen 29 A • D-47559 Kranenburg